Ocular Swan Jacob Autoclavable Gonioprism							
	Product Code	Gonio Mag	Contact OD	Handle Length	Static Field of View		
CE	OSJAG	1.20x	9.5mm	88.17mm	90°	Designed with: Kenneth C. Swan, M.D., University of Oregon, Portland, OR	
C∈	OSJAG8	1.50x	8.0mm	88.17mm	90°		

## **Lens Design**

- The Swan-Jacob Autoclavable Gonioprism is a glass lens designed for easy manipulation for goniotomy and direct viewing gonioscopy for adults and children.
- It provides a clear view of the angle during goniotomy surgical procedures.
- OSJAG8 has 8mm diameter contact suitable for laboratory animals and pediatrics. It has relieved anterior for access to clear cornea.
- It has a black anodized aluminum handle.

## Cleaning & Disinfection

See Cleaning Method 3

## Warnings

- Allow device to air cool approximately 20 minutes before handling and use.
- Use of third party devices and attachments may void the manufacturers' warranty.
- Slight changes to the appearance, such as color, of the adhesive are normal with repeated sterilization and use.

